



STERILIZATION PACKING

HAORUN MEDICAL PRODUCTS CO.,LTD

TEL 0086-574-87706193

WHATSAPP/WECHAT/TEL:008613566542343

EMAIL:SALES@HAORUNMED.COM

WWW.HAORUNMED.COM

STERILIZATION REEL

特卫强 TYVEK



- Tyvek+PET/PE易撕膜
- 高级复合透明薄膜
- 采用环保印刷油墨
- 适用于低温等离子、辐照、环氧乙烷灭菌
- 完全遵照EN 868, ISO11607标准生产
- 适用于医院、实验室、牙科门诊、美容用品等行业的灭菌包装

- Tyvek(1073B/1059B/2FS)
- PET/PE Multilayer Film
- Adopt Environmental Printing Ink
- Suit for Plasma, Gamma and EO
- Comply with EN 868, ISO11607
- Apply to hospital, laboratory, dental clinic, nail&beauty supply, tattoo&piercing supply etc.

Width	
○ 50mm	○ 300mm
○ 75mm	○ 350mm
○ 100mm	○ 400mm
○ 150mm	○ 500mm
○ 200mm	○ 600mm
○ 250mm	
Length	
○ 70m	○ 200m
○ 100m	

平面卷袋 FLAT



- 医用级透析纸(60克或70克)
- 高级复合透明薄膜
- 采用环保印刷油墨
- 带有高压蒸汽、环氧乙烷灭菌指示
- 通过CE, ISO13485认证
- 完全遵照EN 868, ISO11607标准生产
- 适用于医院、实验室、牙科门诊、美容用品等行业的灭菌包装

立体卷袋 GUSSETED



- Medical Grade Paper(60GSM or 70GSM)
- Transparent Composite Film(white, Blue or Green)
- Adopt Environmental Printing Ink
- With STEAM,EO Chemical Indicator
- Approved CE and ISO13485
- Comply with EN 868,ISO11607
- Apply to hospital, laboratory, dental clinic, nail&beauty supply, tattoo&piercing supply etc.

STERILIZATION POUCH

自封灭菌袋 Self Sealing



- 医用级透析纸(60克或70克)
- 高级复合透明薄膜
- 采用环保印刷油墨
- 带有高压蒸汽、环氧乙烷灭菌指示
- 通过CE, ISO13485认证
完全遵照EN 868,
ISO11607标准生产
- 适用于医院、实验室、
牙科门诊、美容用品等
行业的灭菌包装



平面灭菌袋 Heat Sealing/FLAT

Width	Length
○ 50mm	○ 100mm
○ 70mm	○ 130mm
○ 90mm	○ 150mm
○ 100mm	○ 190mm
○ 135mm	○ 230mm
○ 190mm	○ 280mm
○ 350mm	○ 350mm
○ 400mm	○ 500mm
○ 450mm	○ 600mm



立体灭菌袋 Heat Sealing/GUSSETED



- Medical Grade Paper(60GSM or 70GSM)
- Transparent Composite Film(white, Blue or Green)
- Adopt Environmental Printing Ink
- With STEAM, EO Chemical Indicator
- Approved CE and ISO13485
- Comply with EN 868, ISO11607
- Apply to hospital, laboratory, dental clinic,
nail&beauty supply, tattoo&piercing supply etc.

STERILIZATION INDICATOR TAPE/CREPE PAPER/MEDICAL PAPER POUCH

灭菌指示胶带

Sterilization Indicator Tape

- 粘附牢固，易剥离且无残留。
- 适合环氧乙烷、高压蒸汽、等离子灭菌
- 符合标准ISO11607, ISO11140



- Tightly adhering to paper, fabric and plastic
- Easy-removing and no residue
- Suit for EO, STEAM and Plasma Sterilization
- comply with ISO11607 and ISO11140

Width	Length
<input type="radio"/> 12mm	<input type="radio"/> 30m
<input type="radio"/> 19mm	<input type="radio"/> 50m
<input type="radio"/> 25mm	

医用皱纹纸

Medical Crepe Paper



- 采用纯木浆制造（60g或45g）
- 防水，耐撕裂，阻隔细菌
- 可生物降解且无毒无异味
- 符合标准ISO11607
- Made from Virgin Wood Pulp (60GSM or 45GSM)
- Waterproof, tear resistance, bacteria resistance
- Biodegradable, Non-Toxic and Odors
- comply with ISO11607

Width	Length
<input type="radio"/> 30 x 30cm	
<input type="radio"/> 45 x 45cm	
<input type="radio"/> 60 x 60cm	
<input type="radio"/> 75 x 75cm	
<input type="radio"/> 90 x 90cm	
<input type="radio"/> 100 x 100cm	
<input type="radio"/> 120 x 120cm	



- 医用涂胶纸（60或65克）
- 复合透明薄膜（52mic ± 2或医用淋膜纸（58g ± 2）
- 具有良好的阻菌性和透气性，易剥离且热和强度好。
- 采用水溶性环保油墨印刷
- 适合环氧乙烷和伽马射线灭菌
- 适用于一次性医疗器械、无纺布制品和医用敷料等包装。（如输液器、医用纱布、手术衣、乳胶手套等）

- Medical Grade Paper with Coated(60 or 65GSM)
- PE+PET Transparent Film(52mic ± 2) or PE Coated Paper(58GSM ± 2)
- Good bacteria resistance and air permeability, Well sealing strength and excellent peel property.
- Adopt Environmental Printing Ink
- Suit for EO and GAMMA Sterilization.
- Use for packing medical disposable items, non-woven products and rubber products like infusion sets, latex gloves, surgical dressing, gauze etc.

纸纸包装袋 PAPER/PAPER



纸塑包装袋 PAPER/PLASTIC



MEDICAL PACKING MATERIAL

医用涂胶纸

Medical Coated Paper



- 医用原纸涂热熔胶或水性胶
- 附着力强，涂胶均匀
- 良好的阻菌性和透气性
- Coating Hot-Melt Glue or Water-Based Glue
- Well Adhesion and Evenly Coating,
- Good Bacteria Resistance and Air Permeability



Composite Film 复合膜

- 复合牢度强无分层
- 极高的透明度和洁净度
- 种类:
 - 1.三层共挤膜 (PE+PP复合膜)
 - 2.PE+PET复合膜
 - 3.CPP+PET复合膜
- Tightly Compound without Delamination
- Excellent Transparency and Cleanliness
- Type: 1.PE+PP
2.PE+PET 3.CPP+PET

